

PROTESIL ALGICLEAN



Protesil Algiclean is a chemically neutral PH powder which, when mixed with warm water, is excellent in removing deeply engraved and residual alginate from trays and all surfaces in the Dental clinic.

CHARACTERISTICS

- Neutral PH.
- Efficacious even in the most inaccessible areas.
- Ecological formulation, does not contain phosphate or tensio-actives
- Safe to metals due to its neutral pH. Phosphate-free
- Powder concentrated to be mixed with water to obtain 10 lt.
- Peach colour
- Flavour free

PACKAGING

Cod. Prod.	Standard pack
005050	Jar of 500g + measure

VANNINIDENTALINDUSTRIES.r.l.

Via di Campigliano 55/A 50012 GRASSINA (FI) - ITALY Phone +39055644698 Fax +39055 644697

E-mail: info@vanninidental.com www.vanninidental.com