

## **PRESTIGE GYPSCLEAN**



PRESTIGE GYPSCLEAN, Neutral pH cleansing agent for the removal of dental gypsum and gypsum investments from trays, resin prostheses and other instruments.

### TECHNICAL CHARACTERISTICS

- Ready-to-use solution
- Effervescent
- Orange colour
- Action time: 30'-60'

### ADVANTAGES

- Does not contain phosphates or surface-active agents
- Effective even in the least accessible areas
- Gentle

### STORAGE

Keep it securely closed, and store in a cool, dry place at a temperature between 4°C e 28°C (39°F - 82°F). Keep away from heat sources, sparks or flames.

The expiration date printed on the packaging refers to the unopened product still sealed by the production.

The serial number and expiration date are printed on each individual product container. When contacting the distributor, it is important to have this information on hand in order to properly and promptly identify a given product batch.

**IMPORTANT:** Do not use after the expiration date.

### PACKAGING

Product Code/ Description                      058007 - bottle 1 Lt.

VANNINI DENTAL INDUSTRY S.r.l.

Via di Campigliano 55/A 50012 GRASSINA (FI) – ITALY Phone +39 055 644698 Fax +39 055 644697

E-mail : [info@vanninidental.com](mailto:info@vanninidental.com) [www.vanninidental.com](http://www.vanninidental.com)